


## 2.12. Sulforhodamine B (SRB) Assay

CK Cho-Long Kim    HJ Han-Sol Jeong    JM Jung-Soon Mo

Updated date: Apr 14, 2023

 An abbreviated version of this protocol was published in Evidence-based Complementary and Alternative Medicine : eCAM in Feb 2021

Extracts of *Perilla frutescens* var. *Acuta* (Odash.) Kudo Leaves Have Antitumor Effects on Breast Cancer Cells by Suppressing YAP Activity

DOI: 10.1155/2021/5619761

### Detailed protocol

#### SRB staining

##### Materials

- **SRB solution (dye)**: 0.4% SRB (230162, Sigma-Aldrich) in 1% acetic acid diluted with auto ddH<sub>2</sub>O
- **0.1% acetic acid (washing buffer)**: 0.1% acetic acid (695092, Sigma-Aldrich) in auto ddH<sub>2</sub>O
- **10% TCA solution (fixation)**: 10% TCA solution (T0699, Sigma-Aldrich) in auto ddH<sub>2</sub>O
- **10mM Tris base solution (elution)**: Trizma® base (T6066, Sigma-Aldrich) in auto ddH<sub>2</sub>O

##### Procedure

1. Cell seeding in 96well plates (seed a different number of cells depending on their cell sizes)

Example) 1x10<sup>3</sup> cells / each well in 96 well plate

2. Drug treatment: serial dilution (total media has to be 200 ul in each well)

Example) 0.01 – 0.1 – 1 – 2 – 5 – 10 – 20 – 40 – 60 – 80 – 100 – 1000

3. Incubation in 37°C incubator 12 to 16 hours

4. Harvest each plate in once every two days

- 1) Suction media carefully

- 2) Put 100 ul of 10% TCA solution in each well by multichannel pipette

- 3) Incubate the plate in 4°C for 1 h (dark)

- 4) washing TCA solution by D.W about 3 to 5 times

- 5) Add 80ul of 0.4% SRB solution using by multichannel pipette

- 6) Incubate the plate in the dark for 30 min

- 7) Remove the SRB solution with a 0.1% acetic acid buffer (washing)

- 8) Dry for 1h

- 9) Add 150ul of Trizma base solution and shaking 1 h to resolve the staining

- 10) Measure the absorbance in 540 nm

- 11) Analysis the raw data

**How to cite:** (Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Kim, C. , Jeong, H. and Mo, J. (2023). 2.12. Sulforhodamine B (SRB) Assay. Bio-protocol Preprint. [bio-protocol.org/prep2205](https://bio-protocol.org/prep2205).
2. Kim, C., Shin, Y., Choi, S., Oh, S., Kim, K., Jeong, H. and Mo, J.(2021). Extracts of *Perilla frutescens* var. *Acuta* (Odash.) Kudo Leaves Have Antitumor Effects on Breast Cancer Cells by Suppressing YAP Activity. Evidence-based Complementary and Alternative Medicine : eCAM 0(0). DOI: [10.1155/2021/5619761](https://doi.org/10.1155/2021/5619761)

**Copyright:** Content may be subjected to copyright.